

**wapp** Vibe Ball  
Vibration Massage

VERSION: 240103

# Vibe Ball

English\_04

Français\_20

Deutsch\_36

Español\_52

Italiano\_68

# Thank you for purchasing this product

We are sure you will appreciate the care taken in designing and implementing it. We constantly research and combine innovation and technology to create high-quality products.

We apply the most advanced methods and use the materials to always obtain excellent performance. We offer products designed for the most demanding customers who want to achieve excellent results.

**Read these instructions carefully before using the appliance and keep them for future reference.**

## Warnings

**To ensure correct usage of our device, we recommend you read this manual carefully prior to using the product. Please keep the user manual on hand for easy reference.**

**All instructions and precautions outlined in this manual should always be followed when using the product.**

1. This product contains batteries that are not replaceable by the user.
2. Batteries are only replaceable by skilled persons.
3. Do not leave the device unattended when connected to mains power.
4. When not in use and before removing or replacing parts, the device must be unplugged from the mains power outlet.
5. Do not operate the device under a pillow or blanket/bedding, as overheating may occur causing a fire, electric shock or injury to persons.
6. Do not carry this device by its supply cord or use the cord as a handle.
7. Do not operate the device whilst blocking the air openings. The air openings must remain dust, lint and hair free.
8. Do not use the device in damp rooms such as bathrooms or laundries, when showering or bathing.
9. Supervision is required if the device is being operated by or near children, or disabled persons.
10. Do not place the device's cord or the device itself, near heated surfaces.

## Warnings

11. Do not use the device if it has been dropped, immersed in water or appears damp.
12. The device is not waterproof; do not expose it to rain or moisture
13. Do not use on areas with swelling, burns, inflammation, open wounds, fresh wounds, hematomas, bleeding or other injuries
14. Do not use the product if it is not working properly, dropped or damaged
15. If you do feel uncomfortable while using it, you must try to change your position or contact pressure. If you still feel uncomfortable, stop using it.
16. The appliance is designed for private use. It's not suitable for commercial use, e.g., in hairdressing salons, gyms, medical purposes or beauty studios.
17. This device is not intended for use by persons (including children) with limited physical, sensory or mental abilities, or persons lacking experience/or knowledge, unless any such person is supervised by a person responsible for his or her safety.

## Warnings

This symbol indicates that this product should not be treated as domestic waste. Once its useful life has ended, it must be taken to a relevant collection point for the recycling of electrical appliances. Through the correct recycling of batteries and electrical devices, we contribute to avoiding risks to environmental health and safety.

**Product Disposal Note:** If this product is damaged beyond repair or, for any reason, its useful life is deemed to have expired, please inform yourself about local, state and European regulations regarding the proper disposal and recycling of products containing lead, batteries, plastics, among other materials, as well as collection points for these types of products.

## Warnings

Li-ion

These symbols indicate that this product contains batteries and cannot be disposed as domestic waste. Once its useful life has ended, it must be taken to a relevant collection point for the correct recycling of its batteries.

This product complies with the CE Marking. You can find the European Declaration of Conformity at the web link provided: **<https://www.weelko.com>**

The digital manual for this product is available at the web link provided: **<https://www.weelko.com>**

## **Cleaning / Care**

- Turn the device off after using and before cleaning
- Clean the device with a damp cloth only
- Liquid should not enter the device in any circumstances
- Do not use the device until it is completely dry after cleaning
- Do not clean the device in a dishwasher
- Do not use abrasive cleaning products or hard brushes
- If you don't plan to use the device for a long period of time, store it in the original packaging in a dry environment

## **Before using the device**

- Do not use the device if the packaging, the device itself, or the USB power cord show signs of damage.
- The device is not to be used for commercial use or medical purposes.
- Prior to using the device, please consult your medical health practitioner if you have any health concerns.

## Package content


- 1 Vibration massage Ball
- 1 Type C charging cable
- 1 user manual

## Features


- Heat therapy
- 360° massage ball
- 4 vibration modes
- Cordless design & lightweight
- Deep tissue massage
- One-hand operation
- High-intensity vibration
- Auto shut down
- Full body massage
- Built-in rechargeable battery

## Specifications

- 11 cm diameter
- Battery 1.200mAh
- Input: DC 5V/1A
- 4 vibration modes: Low-Medium-High-Auto
- Material: Silicon Rubber & ABS
- 7 indicators: 3 for vibration modes, 3 for battery and 1 for heating
- Charge time: 2,5 hours
- Using time: Around 60-120 minutes depending on the selected vibration mode and if heating mode is working.

## Parts of the product

1. Charging point
2. Battery indicator
3. Heat function indicator
4. Center button
5. Vibration mode indicator

## Battery charging

The device must be fully charged prior the first use.

It takes approximately 2,5h to fully charge. The charging cable needs to be connected to a suitable USB output (5V/1A).

While charging, if the battery charge is below 30%, one light will flash. When the battery is charged around 30 - 70%, two lights will flash. When the battery is charged around 70 - 99%, three lights will flash. When fully charged, 3 lights will illuminate without flashing. Then the device can be used for approximately 1 to 2 hours, depending on speed level used and heating.

## Quickstart guide

Once the device is ON, the battery indicator will show the level of remaining battery.

- 3 LEDs ON: 70%–100% of battery
- 2 LEDs ON: 30%–70% of battery
- 1 LED ON: under 30% of battery

**The vibration mode indicator will illuminate 1, 2 or 3 lights depending on the mode selected:**

- Press the center button to turn ON the device. Low-intensity massage mode will start working, the first blue indicator will turn on.
- Press the center button a second time to change to medium-intensity massage, the second blue indicator will turn on.
- Press it a third time to switch to high-intensity massage, the third blue indicator will turn on.
- Press it fourth time to switch to automatic intensity massage, three blue indicators will turn on.
- When you press the button for a fifth time, the device will turn OFF.

## Heating function

- Long press the center button for three seconds to turn ON/OFF the heating function.
- This function can be activated at any vibration model.
- The heating function will remain active until it's deactivated or until the device is turned OFF.

## Safety feature

- The device has a safety feature which automatically turns OFF after 15 minutes of use.
- It can be turned back ON again and used until recharging is required.

## Note

- The vibrating ball is designed with foam ridges around the surface. The ridges help to relieve muscle tension and relaxes the areas you are targeting. The more pressure you put on the ball the deeper the massage will be.
- Slowly roll your body back and forth along the vibration device. Use controlled movements, and pay attention to how your muscles respond to the vibrations. Roll each body area for 1-2 minutes.
- It is recommended to start using the lowest speed and increase gradually the vibration mode depending on sensations during the massage.
- If you experience intense pain, stop applying pressure to the device or turn OFF the vibration function.

## Applications and benefits

Use the device by rolling it over the muscles as indicated in the images below.

Hips

Waist


Arms


Legs

Back

Calves

### Benefits:

The device may assist:

- Release tight muscles and tension
- Boost circulation
- Reduce stress
- Improves relaxation and sleep

## Troubleshooting

| PROBLEMS                                   | CAUSE                                                                                               | SOLUTION                                                                                                                                       |
|--------------------------------------------|-----------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Device does not work:</b>               | The device is either not charged up sufficiently, or the charger plug may not be inserted properly. | Charge the device with the USB cable supplied.                                                                                                 |
|                                            | Device not switched on.                                                                             | Press the On button.                                                                                                                           |
|                                            | The device will automatically switch off after 15 minutes of use.                                   | Press the On button again.                                                                                                                     |
| <b>Device gets hot, or cannot be used:</b> | The product has been used repeatedly.                                                               | If the device is being charged, disconnect the charger and allow it to cool. Always disconnect the device from the charger when fully charged. |
|                                            | If the 15 minutes continuous use time has been exceeded, the device can overheat.                   | Turn the device off and allow to cool down for 30 minutes.                                                                                     |

# Warranty


**Thank you for purchasing your new Wapp Vibe Ball body massager from Weelko.**

Weelko guarantees that this product will be free from defects in materials and workmanship for 12 months from the date of purchase, provided that the product is used in accordance with the accompanying recommendations or instructions.

Please keep your receipt as proof of purchase and contact the Weelko Support Team at **info@weelko.com** or by calling **+34 938249020**. Warranty claims and expenses incurred in returning this product can be directed to the Weelko Support Team. **Barcelona SL, Camí de Garreta 92 · 08600 Berga, Spain.**

**Acquisition date:**

**Signature and seal of the sales establishment:**

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