

# Directa Crowns Molars Mould Guide

Translucent/Translucide/Translucida/Translucenti

## Upper Right/Höger ÖK/Haut droit/OK-Rechts/Superior Derecha/Superiore Destro

mm ↓	9.2	9.5	9.6	10.4	11.0	11.4	12.3	12.9	13.3
mm ↔	9.9	10.3	10.7	11.1	11.5	11.9	12.6	13.2	13.8



## Upper Left/Vänster ÖK/Haut gauche/OK-Links/Superior Izquierda/Superiore Sinistro

mm ↓	10.0	10.2	10.8	11.3	11.7	12.2	12.6	13.4	14.0
mm ↔	9.9	10.3	10.7	11.1	11.5	11.9	12.6	13.2	13.8



## Lower Right/Höger UK/Bas droite/UK-Rechts/Inferior Derecha/Inferiore Destro

mm ↓	7.9	8.2	8.8	9.2	9.5	10.0	10.7	11.1	11.6
mm ↔	9.9	10.3	10.7	11.1	11.5	11.9	12.6	13.2	13.8



## Lower Left/Vänster UK/Bas gauche/UK-Links/Inferior Izquierda/Inferio Sinistro

mm ↓	8.1	8.3	8.9	9.1	9.2	9.5	10.2	10.5	11.1
mm ↔	9.9	10.3	10.7	11.1	11.5	11.9	12.6	13.2	13.8



**DIRECTA**  
directadental.com



# Directa Crowns Mould Guide

Opaque/Opak/Opaca/Opache

Translucent/Translucide/Translucida/Traslucenti

Centrals/Centraler/Centrales/OK-Front/Centrales/Centrali															
mm	↓	14.0	13.0	12.4	12.0	11.5	11.0	10.6	10.6	11.0	11.5	12.0	12.4	13.0	14.0
mm	↔	10.0	9.2	8.7	8.4	8.2	7.8	7.5	7.5	7.8	8.2	8.4	8.7	9.2	10.0
<p>100 101 10 11 12 13 14 15 16 17 18 19 102 103</p> <p>centrals</p> <p><b>DIRECTA</b> selecta</p>															

  

Laterals/Lateraler/Latérales/OK-Front/Laterales/Laterali															
mm	↓	11.0	10.4	10.1	9.5	8.8	8.2	8.2	8.8	9.5	10.1	10.4	11.0		
mm	↔	7.3	6.8	6.5	6.2	5.9	5.5	5.5	5.9	6.2	6.5	6.8	7.3		
<p>2 20 21 22 23 24 25 26 27 28 29 200</p> <p>laterals</p> <p><b>DIRECTA</b> selecta</p>															

  

Cuspids/Hörntänder/Canines/Eckzähne/Caninos/Cuspidi															
mm	↓	13.0	12.4	11.6	11.4	11.1	10.8	10.5	10.5	10.8	11.1	11.4	11.6	12.4	13.0
mm	↔	9.0	8.5	8.1	7.9	7.7	7.5	7.3	7.3	7.5	7.7	7.9	8.1	8.5	9.0
<p>300 301 30 31 32 33 34 35 36 37 38 39 302 303</p> <p>cuspids</p> <p><b>DIRECTA</b> selecta</p>															

  

Upper bicuspids/ÖK-premolaren/Prémolaires supérieures/OK Prämolaren/Premolares superiores/Bicuspidi superiores															
mm	↓	9.8	9.4	8.9	8.6	8.4	7.5	8.1	8.5	9.0	9.5				
mm	↔	7.4	6.9	6.7	6.3	6.0	6.0	6.3	6.7	7.2	7.4				
<p>40 41 42 43 44 45 46 47 48 49</p> <p>bicuspids upper</p> <p><b>DIRECTA</b> selecta</p>															

  

Lower bicuspids/UK-premolaren/Prémolaires inférieures/UK Prämolaren/Premolares inferiores/Bicuspidi inferiores															
mm	↓	9.8	9.5	9.2	8.8	8.2	7.1	7.3	8.0	8.1	8.4				
mm	↔	7.4	7.2	6.8	6.6	6.3	6.4	6.8	7.1	7.5	8.0				
<p>50 51 52 53 54 55 56 57 58 59</p> <p>bicuspids lower</p> <p><b>DIRECTA</b> selecta</p>															

  

Lower incisors/UK-incisiver/Incisives inférieures/UK Front/Incisivos inferiores/Incisivi inferiori															
mm	↓	12.0	11.5	11.0	10.5	10.0	8.0	8.5	9.1	9.6	10.3				
mm	↔	6.2	5.6	5.3	5.1	4.9	4.5	4.9	5.2	5.5	5.8				
<p>60 61 62 63 64 65 66 67 68 69</p> <p>incisors lower</p> <p><b>DIRECTA</b> selecta</p>															

